

Blades Biological Ltd

COWDEN . EDENBRIDGE . KENT . U.K. . TN8 7DX

Telephone: (01342) 850242 Fax: (01342) 850924

Email: info@blades-bio.co.uk Online: <http://www.blades-bio.co.uk/>

Horse Blood

Oxilated and Defibrinated are both treated to prevent clotting.

Oxilated contains potassium oxilate which dilutes calcium and prevents it taking its place in the clotting process.

Defibrinated blood has had its fibrin removed mechanically making it unavailable to form a clot.

Defibrinated blood is the best for centrifuging, with the cell forming on the bottom and the plasma on the top.