

NATURES NUTRITION

25101 COUNTY ROAD 35
WRAY, CO. 80758

DATE: July 22, 2009
ITEM: A14000F

CERTIFICATE OF ANALYSIS

PRODUCT NAME: CAROLINA BIOLOGICAL ALFALFA TABLETS

LOT NUMBER: G170209 P.O. NUMBER: 509062 ITEM CODE: 14578

I. CONTENTS:

Alfalfa Leaf	654.1 mg = 10 grains
Calcium Sulfate	70.6 mg
Silica	12.8 mg
Croscarmellose Sodium	7.5 mg
Total Weight	745.0 mg

II. PHYSICAL:

Weight	745 mg ± 5%
Diameter	1/2 inch round
Thickness	0.242 ± .010 inches
Hardness	8 ± 2 kg/cm
Friability	1% maximum
Disintegration	30 minutes maximum

III. MICROBIOLOGICAL RESULTS:

Aerobic Plate Count	<10,000 cfu/g
Yeast & Mold	< 100 cfu/g
E. coli	Negative
Salmonella spp.	Negative
Staphylococcus aureus	Negative

ALVA DETERDING
TECHINCAL DIRECTOR